
SkinElly

MESO

Hyaluronic acid gel

INSTRUCTION FOR MEDICAL USE OF MEDICAL PRODUCT

Medical Product Name:

Intradermal Implant «SkinElly» with Sodium Hyaluronate, SkinElly MESO Version 24 mg/ml composed of:

- 1 ml syringe with an intradermal implant with sodium hyaluronate - 1 pc.
- Needles - 0,3 x 13 mm (30G) - 2 pcs.
- Self-adhesive labels (sticker) - 3 pcs.
- Instruction for use (package leaflet) – 1 pc.

MANUFACTURER/DEVELOPER OF THE MEDICAL PRODUCT:

ACTIVE MEDICAL PRODUCTS LIMITED, Suite 1,
5 Percy Street, Fitzrovia, London W1T 1DG (UK)

CORRECT USE OF THE MEDICAL PRODUCT:

Intradermal Implant «SkinElly» with Sodium Hyaluronate, SkinElly MESO Version 24 mg/ml is indicated for use in cosmetology, plastic and aesthetic surgery for the correction of involuntal skin changes and non-surgical intimate contour plastic surgery.

Intradermal Implant «SkinElly» with Sodium Hyaluronate, SkinElly MESO Version 24 mg/ml is an injection, transparent, colorless, sterile, pyrogen-free, fully biosynthetic, 100% bioabsorbable in human body tissues, placed in a sterile glass syringe with a volume of 1.0 ml pre-filled with a colorless, transparent and viscoelastic substance.

INGREDIENTS:

Hyaluronate sodium 24 mg/ml (2,4 %)

Phosphate buffer solution

INDICATIONS FIELD:

Intradermal Implant «SkinElly» with Sodium Hyaluronate, SkinElly MESO Version 24 mg/ml is used in plastic and aesthetic surgery, aesthetic medicine and cosmetology for the correction of involuntal skin changes. It is to be injected into the superficial and mid dermis to correct small and medium wrinkles, skin defects associated with dehydration, decrease in turgor, skin tone and elasticity. It increases the regenerative potential of the skin, due to the high degree of hygroscopicity provides a powerful and sustained moisturizing effect, improves the skin structure and color. Provides an expressed rejuvenating and remodeling effect. It is a perfect solution for patients with a rugulose type of aging. Prolonged hydration of the skin of the face, neck, décolleté and back of the hands. Atrophic scars (post-acne).

It is used for the non-surgical intimate contour plastics: xerosis of the anogenital zone, sclerosing lichen, leukoplakia, neoplasias of the cervix and vulva, dyspareunia, vulvodynia, vaginismus, fissures, chronic pelvic pain, restoration of natural moisturizing of the vaginal mucosa and anogenital zone, with complex treatment of herpes.

CONTRAINDICATIONS TO THE MEDICAL PRODUCT USE:

General:

- Medical products only for intradermal use;
- Allergic reaction to any components that make up the medical product;
- Pregnancy and lactation;
- In combination with laser procedures or peeling;
- Blood vessels administration;
- Systemic treatment with retinoids;
- Systemic anticoagulant and antithrombotic therapy;
- Systemic immunosuppressive therapy;
- Inflammatory skin infiltration;
- Progressive chronic dermatosis (Koebner Phenomenon);
- Infection caused by the herpes virus in the acute stage;
- Paraneoplastic process at the administration site;
- Severe somatic disease at the decompensation stage;
- Acute infectious diseases, feverish conditions;
- Acute exacerbation of chronic diseases, accompanied by decompensation processes of the II-III stage;
- Autoimmune diseases;
- Long-term treatment at the dentist (dental prosthetics) during the procedure;
- Under the age of 18 (without the consent of the parents or their substitute);
- Hypersensitivity to hyaluronic acid drugs in past medical history;
- Bleeding disorder: hemophilia;
- Porphyria.

Absolute:

- Tendency to form hypertrophic and keloid scars or any other impaired healing;
- Permanent fillers in the administration area.

Relative:

- Inflammation and infectious processes in the intended administration area (acne, herpes, eczema, etc.);
- Violations of the epidermis integrity in the administration zone (for example, post-pilling period, etc.).

CONTRAINDICATIONS FOR USE IN NON-SURGICAL INTIMATE CONTOUR PLASTIC SURGERY:

- Hypersensitivity to any component of a medical product;
- Burdened allergic history;
- Diabetes mellitus in decompensated stage;
- Tuberculosis;
- Clinically significant bleeding disorders or anticoagulants use;
- Malignant neoplasms;
- Pregnancy or breastfeeding;
- Chronic diseases at the decompensation stage;
- The active phase of herpes infection;
- Infectious diseases of the urogenital tract in the acute phase;
- Mental illness;
- Urinary bladder emptying disorder;
- Acute infectious and inflammatory diseases of the urogenital and extragenital systems and exacerbations of chronic diseases;
- Period of menstruation;
- Infectious-allergic bronchial asthma;
- Presence of permanent fillers in the correcting zones;
- Autoimmune diseases;
- Tendency to develop hypertrophic scars, acute inflammatory processes at the intended injection sites;
- Inflammatory processes and sexually transmitted diseases;
- Individual intolerance to hyaluronic acid;
- Presence of an unknown filler or drugs of other physical and chemical properties of another manufacturer in the planned projection zone that was injected less than 12 months ago;
- Age up to 18 years.

PRECAUTIONS:

There were no studies on the use of a medical product depending on the form of execution of the Intradermal implant with «SkinElly» sodium hyaluronate in patients with already implanted medical products based on hyaluronic acid. It is necessary to determine the "benefit – risk" ratio of the proposed procedure before using medical product.

- Check the expiration date and integrity of the medical product packaging before use;
- Use only for one patient;
- Do not re-sterilize;
- Comply with sanitary and epidemiological requirements during injection procedures as there is a risk of infection;

- Do not expose to intense heat (sunshine or tanning beds) or significant cooling of the administration site;
- Dispose of the syringe with implant residues after use.

Note: Removing or replacing an **intra**dermal implant is not possible. The use of various hardware techniques (for example, d'Arsonval, electroporation, microcurrent therapy, vaporization) or the administration of the enzyme hyaluronidase or Longidase is possible for removal or accelerated absorption at the discretion of the practitioner.

INTERACTIONS WITH OTHER MEDICAL PRODUCTS AND MEDICINAL DRUGS

The use of a medical product Intradermal implant «SkinElly» with sodium hyaluronate in combination with other medical products and medicinal drugs has not been studied.

Note 1: Removing or replacing an Intradermal implant «SkinElly» with sodium hyaluronate is not possible. The use of various hardware techniques (for example, d'Arsonval, electroporation, microcurrent therapy, vaporization) or the administration of the enzyme hyaluronidase or longidase is possible for removal or accelerated absorption at the discretion of the practitioner.

Note 2: The effect of the Intradermal implant «SkinElly» with sodium hyaluronate on vision and reaction time has not been analyzed and therefore there are no special instructions for drivers and/or persons operating complex technological devices.

The doctor should not:

- Use together with other drugs, both in one syringe, and in one procedure;
- Administer to patients with known bleeding disorders or patients receiving medicinal drugs that increase blood clotting time (vitamin K antagonists, heparin, or immunosuppressant drugs), as they may increase the risk of bleeding at the site of administration;
- With unsatisfactory experience of this type of injection;
- Use before or after laser treatment or chemical peeling, as this may cause inflammation of the injected area.

The patient should not:

- Subject the treated area to heat (sunbathing, sauna, steam bath, etc.) or extreme cold until any signs of local inflammation disappear.

There is a risk of infection characteristic of all subcutaneous injections. Therefore, standard precautions are used for this type of procedure when administering the product.

You must ensure that the individual sterile packaging is not damaged, the expiration date of the syringe and needles has not expired before use.

Do not use the product if the contents of the syringe look cloudy.

It is forbidden to administer into blood vessels, bones, ligaments, tendons and muscles and into birthmarks.

The remains of the medical product must be disposed of!

It is forbidden to use medical product:

- With sterile packaging integrity breach. Do not re-sterilize!
- If the syringe cap of the Luer-Lok cap is gapping;
- Expired;
- With color changes or transparency of the syringe contents.

EXPECTED AND PREDICTABLE ADVERSE REACTIONS

Transient reactions:

As a manifestation of transient reactions, the appearance of short-term bleeding (1-2 minutes) and the appearance of bruising, pain and itching at the moment of injection for a short time (10-15 minutes) are possible. Perhaps the appearance of erythema, as a result of local irritating or compensatory action, which, in turn, passes quite quickly (no more than 2 days). The manifestation of the injection site mass (for no more than 4-5 days), discoloration of the skin and an increase in sensitivity at the injection site

Preservation of any reactions over a specified time and (or) the addition of inflammatory symptoms (temperature, infiltration, fluctuation, puffiness and swelling, granuloma, etc.) indicate complications and require that the patient be treated by a specialist with appropriate profiles using pharmaceuticals drugs.

BIODEGRADATION

Biodegradation of the Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly MESO 24 mg/ml depends on the design (concentration) and the individual characteristics of the patient:

BIODEGRADATION OPTIONS

Hydrolytic degradation: The destruction of the material in the aquatic environment under the influence of:

- acids and alkalis (non-enzymatic hydrolysis);
- enzymes (enzymatic hydrolysis);
- reactive oxygen intermediates and free radicals (oxidative destruction, which is significantly accelerated in the presence of metal ions as a catalyst).

Cell degradation: Biomaterial destruction by macrophages and mast cells.

Bacterial degradation: The presence of bacteria in the area of implant insertion can accelerate biodegradation due to enzymatic and oxidative degradation.

Mechanical degradation: Fragmentation and/or change of the implant shape under the action of static or dynamic loads can accelerate the destruction of the material (active facial expressions, massage, strokes, etc.).

ELIMINATION PATHWAYS:

Sodium hyaluronate is eliminated from the body depending on the individual characteristics of the patient; in the human body sodium hyaluronate is laminated by hyaluronidase enzymes. The decomposition products of sodium hyaluronate are oligosaccharides and extremely low molecular hyaluronates, which are naturally eliminated from the human body through the liver and kidneys.

THE TIME OF ONSET OF THERAPEUTIC EFFECT

The time of onset of the therapeutic/cosmetological effect depends on the individual characteristics of the patient.

THE TIME OF THERAPEUTIC EFFECT PRESERVATION

The time of therapeutic effect preservation depends on the individual characteristics of the patient and in the human body.

INFORMATION ON THE PROCEDURE OF WORKING WITH A MEDICAL PRODUCT (PROCEDURE AND/OR USAGE RULES)

Preparation for Use

The procedure is carried out in accordance with the standard medical practice and local regulations. The product should be used by doctors who have received special training in injection techniques for filling skin folds for a successful procedure.

Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly MESO Version 24 mg/ml to be used as it is delivered.

Carefully inspect all parts for damage.

- Before injection, carefully inspect the SkinElly implant visually. You may need to perform several additional procedures to achieve optimal results. In order to enable effective evaluation of the implantation results, a break should be made between the courses of treatment for at least seven days.
- Before starting treatment, patients should be informed about the indications for the product use, contraindications, incompatibility and potential side effects. Before the procedure, you should read the full history of the patient and carefully examine the area to be treated.
- You should also evaluate the patient's ability to cope with pain and, if necessary, apply the most appropriate form of anesthesia. In the case of administration of an anesthetic substance, ice should be applied to the treated area to reduce swelling and local stretching.

- Before starting the procedure, thoroughly clean the area to be treated with soap and water and disinfect with an alcohol swab.

Intradermal Implant «SkinElly» with Sodium Hyaluronate, SkinElly MESO Version 24 mg/ml: is to be injected into the superficial and mid dermis

Injection zones:

Face, neck, décolleté, hands, earlobes, limb skin, trunk skin

Injection technique:

- Micropapules.
- Microdroplet to the area of «kiset» wrinkles, «proud» wrinkles.
- Microbolus in the red border of lips, earlobes.
- Linear-retrograde with separation of superficial linear wrinkles.

Injection zones and injection techniques for non-invasive intimate contour plastic surgery:

Superficial and deep layers of the dermis of the anogenital region or the intramucosal layer of the vagina, vulva.

Method of injection:

Papular and micropapular injection technique.

Immediately stop the procedure if you suspect a puncture of the vessel.

If significant resistance is felt during the gel injection, the needle should be moved deeper to ensure easier insertion of the implant. If resistance remains, try inserting the implant from another entry point. In case of continuous resistance, it may be necessary to replace the injection needle and/or syringe.

Any excessive pressure on the syringe should not be applied at any time during the procedure.

If blanching occurs immediately after the gel injection, discontinue the injections and massage the treated area until the skin returns to its normal color. If the skin has not regained its normal color, the product administration process should not be resumed.

Superficial injections or the injection of large doses of the product may cause discoloration and nodulation.

The duration of the correction of the Intradermal implant with sodium hyaluronate «SkinElly» depends on the choice of performance, the area of administration and the characteristics of the patient. During the procedure after the injection, a light massage is desirable for better redistribution in the tissues.

INSTRUCTIONS FOR THE PATIENT

- Patients should adhere to the following recommendations: avoid excessive physical exertion, as well as exposure to direct sunlight and tanning lamps or extreme weather conditions for 24 hours.

- Patients should be informed of the application of a wrap with ice or a cold compress to the treated area within 24 hours after the procedure.
- Patients are also encouraged to talk less, smile and laugh for one week after the procedure.
- Patients who underwent oral procedures should rinse mouth with saline every 3–4 hours, within one week after the procedure.
- If necessary, patients can use oral analgesics as recommended by their physician.
- If the inflammatory process lasts more than a week or other side effects appear, you should immediately consult a doctor.
- You should refrain from applying makeup for 12 hours after the injection, you should avoid prolonged contact with sunlight, ultraviolet light, gels, do not go to the sauna or the bath for 2 weeks after the injection.
- Do not consume alcoholic beverages, anti-inflammatory drugs, anticoagulants and drugs accelerating the metabolism before and after the procedure without consulting a doctor.
- Patients who have undergone non-surgical intimate body sculpting are advised to avoid any mechanical, chemical or thermal effects in the area of the procedure for the next 14 days.

THE OPERATING CONDITIONS OF THE MEDICAL PRODUCT:

The Intradermal implant «SkinElly» with sodium hyaluronate can only be used by medical personnel with appropriate qualifications and familiar with injection techniques.

- Temperature from + 32 °C to +42 °C;
- Humidity from 10 to 80% (non-condensing);
- Atmospheric pressure from 500 to 1060 hPa.

Avoid direct sunlight, freezing, and condensation.

INCOMPATIBILITY

Since hyaluronic acid can cause the release of a disinfectant tertiary (medium) ammonium, for example, benzalkonium chloride or chlorhexidine, the drug should not be near these substances or come into contact with the devices with which these substances were used.

Upkeep and periodic maintenance, including cleaning and disinfection

In the production process of SkinElly, hyaluronic acid undergoes several steps of purification from impurities of residual proteins and endotoxins peculiar to bacterial cell.

Due to the low concentration of residual ingredients of bacterial origin, SkinElly does not exhibit allergenic properties and allergy test before the procedure is not required.

STORAGE CONDITIONS AND TRANSPORTATION STORAGE AND ACTIONS AFTER USE

After completion of treatment, the syringe must be disposed of in accordance with established rules. You must dispose of the syringe with the remnants of the gel and the needle, as they are potentially infected.

STORAGE CONDITIONS

For Intradermal Implant «SkinElly» with Sodium Hyaluronate, SkinElly MESO Version 24 mg/ml

- Temperature from + 2 to +20 °C
- Relative humidity 35-85%;
- Atmospheric pressure 500 - 1060 hPa.
- Non-condensing, away from sunlight.

Ensure packaging integrity before use.

TRANSPORTATION CONDITIONS

Transport at + 2 to +20 °C, away from sunlight, in hermetically sealed packaging.

DISPOSAL OF A MEDICAL PRODUCT

Disposal of a medical product Intradermal implant «SkinElly» with sodium hyaluronate, SkinElly MESO version 24 mg/ml in the territory of the Russian Federation is carried out according to SanPiN 2.1.3684-21 belongs to Class B (epidemiologically hazardous waste) and is a potentially infectious waste.

The used medical product, including needles for injection, not fully used syringes are disposed of according to the specified hazard class of medical waste according to SanPiN 2.1.3684-21. Expired medical products, medical products in damaged packaging should be disposed as class A of medical waste according to SanPiN 2.1.3684-21.

THE SHELF LIFE OF THE MEDICAL PRODUCT

The shelf-life is 2 years from the moment of sterilization while meeting the conditions of transportation and storage.


The opened medical product should be used immediately and disposed of immediately after use.

OFFICIAL REPRESENTATIVE CONCERNING THE MEDICAL PRODUCTS IN THE RUSSIAN FEDERATION

ACTIVE MEDICAL PRODUCTS Limited Liability Company

Russia, 123100, Moscow, 16 Shmitovsky Passage, bld. 2,
Room 1, Office 13, Floor 2

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